

STAT6 Rabbit mAb


ABZ2501

Product Overview

Name	STAT6 Rabbit mAb
Accession(Primary)	P42226
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 94 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; IP-1:20-1:50

Target information (P42226)

Gene Name	STAT6
About Gene Symbol	STAT6 (Signal Transducer And Activator Of Transcription 6) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of STAT6 in Jurkat, CHO-K1, HeLa lysates using STAT6 antibody.