

STAT1 Rabbit mAb


ABZ2502

Product Overview


Name	STAT1 Rabbit mAb
Accession(Primary)	P42224
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:100-1:1000; IP-1:10-1:100

Target information (P42224)


Gene Name	STAT1
About Gene Symbol	STAT1 (Signal Transducer And Activator Of Transcription 1) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cell



Western blot analysis of STAT1 in MCF-7 lysates using STAT1 antibody.



Western blot analysis of STAT1 in MCF-7, HepG2, U2OS, U251MG, K562 lysates using STAT1 antibody.



Immunocytochemistry analysis of STAT1(green) in Hela using STAT1 antibody, and DAPI(blue)