

STAT1 Rabbit mAb


ABZ2503

Product Overview


Name	STAT1 Rabbit mAb
Accession(Primary)	P42224
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Recommended Dilution	WB-1:4000-1:40000; IP-1:20-1:50

Target information (P42224)

Gene Name	STAT1
About Gene Symbol	STAT1 (Signal Transducer And Activator Of Transcription 1) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cell



Western blot analysis of STAT1 in Raw264.7, mouse lung lysates using STAT1 antibody.



Western blot analysis of STAT1 alpha in K562, rat Brain, Hela lysates using STAT1 alpha antibody.