

CD86 Rabbit mAb


ABZ2504

Product Overview


Name	CD86 Rabbit mAb
Accession(Primary)	P42081
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:2000-1:20000; IHC-1:100-1:500; ICC/IF-1:200-1:500; IP-1:20-1:50

Target information (P42081)


Gene Name	CD86
About Gene Symbol	




Western blot analysis of CD86 in Daudi/C6 lysates using CD86 antibody.



Western blot analysis of CD86 in Raji lysates using CD86 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CD86 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of CD86 in K562 using CD86 antibody.