

ACTR1B Rabbit mAb


ABZ2505

Product Overview


Name	ACTR1B Rabbit mAb
Accession(Primary)	P42025
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:200-1:500; IP-1:50-1:100

Target information (P42025)


Gene Name	ACTR1B
About Gene Symbol	



Western blot analysis of ACTR1B in HeLa, A549, HL-60, U251, U87-MG lysates using ACTR1B antibody.



Immunocytochemistry analysis of ACTR1B (green) in HeLa using ACTR1B antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using ACTR1B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.