

CD63 Rabbit mAb

ABZ2506


Product Overview

Name	CD63 Rabbit mAb
Accession(Primary)	P41731
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 26 kDa; Observed MW: 40-70 kDa
Recommended Dilution	WB-1:2000-1:20000; IHC-1:200-1:400; IF: 1:50


Target information (P41731)

Gene Name Cd63

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded mouse kidney using CD63 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD63 in Raw264.7
lysates using CD63 antibody.