

4-1BBL Rabbit mAb

ABZ2507


Product Overview

Name	4-1BBL Rabbit mAb
Accession(Primary)	P41273
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 27 kDa; Observed MW: 25-32 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P41273)


Gene Name TNFSF9

About Gene Symbol



Western blot analysis of 4 1BBL in 293T,Raji lysates using 4 1BBL antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of TNFSF9 in HEK293
lysates using 4 1BBL antibody.