

MMP9 Rabbit mAb

ABZ2508

Product Overview

Name	MMP9 Rabbit mAb
Accession(Primary)	P41245
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 81 kDa; Observed MW: 92 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20

Target information (P41245)

Gene Name Mmp9**About Gene Symbol**

