

HNF 4 alpha Rabbit mAb


ABZ2509

Product Overview

Name	HNF 4 alpha Rabbit mAb
Accession(Primary)	P41235
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (P41235)

Gene Name	HNF4A
About Gene Symbol	



Western blot analysis of HNF4alpha in A549, MCF-7 lysates using HNF 4 alpha antibody.