

Prostaglandin D Synthase Rabbit mAb


ABZ2510

Product Overview



Name	Prostaglandin D Synthase Rabbit mAb
Accession(Primary)	P41222
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 21 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (P41222)


Gene Name	PTGDS
About Gene Symbol	



Immunocytochemistry analysis of Prostaglandin D Synthase (green) in K562 using Prostaglandin D Synthase antibody, and DAPI (blue).



Western blot analysis of Prostaglandin D Synthase in K562, C6, 3T3, HeLa lysates using Prostaglandin D Synthase antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using Prostaglandin D Synthase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.