

Aquaporin 2 Rabbit mAb


ABZ2512

Product Overview

Name	Aquaporin 2 Rabbit mAb
Accession(Primary)	P41181
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; IP-1:20-1:50

Target information (P41181)

Gene Name	AQP2
About Gene Symbol	



Western blot analysis of Aquaporin 2 in Raoms,
Raji lysates using Aquaporin 2 antibody.