

STAT3 Rabbit mAb


ABZ2517

Product Overview


Name	STAT3 Rabbit mAb
Accession(Primary)	P40763
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Hamster
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:20-1:50; FC-1:20; IP-1:20

Target information (P40763)

Gene Name	STAT3
About Gene Symbol	STAT3 (Signal Transducer And Activator Of Transcription 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of STAT3 in Jurkat, C6, CHO-K1 lysates using STAT3 antibody.



Immunocytochemistry analysis of STAT3(green) in Hela using STAT3 antibody, and DAPI(blue)