

MLH1 Rabbit mAb


ABZ2518

Product Overview


Name	MLH1 Rabbit mAb
Accession(Primary)	P40692
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:20-1:50; IF-1:20-1:50; IF-1:20; IP-1:20

Target information (P40692)

Gene Name	MLH1
About Gene Symbol	



Western blot analysis of MLH1 in K562, C6, 3T3, HeLa lysates using MLH1 antibody.



Immunocytochemistry analysis of MLH1 (green) in HeLa using MLH1 antibody, and DAPI (blue)