

Thrombopoietin Receptor Rabbit mAb

ABZ2520


Product Overview

Name	Thrombopoietin Receptor Rabbit mAb
Accession(Primary)	P40238
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P40238)

Gene Name MPL

About Gene Symbol



Western blot analysis of TPOR/MPL in K562, C6 lysates using Thrombopoietin Receptor antibody.

