

MMP12 Rabbit mAb

ABZ2521


Product Overview


Name	MMP12 Rabbit mAb
Accession(Primary)	P39900
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54,45 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (P39900)

Gene Name MMP12

About Gene Symbol





Western blot analysis of MMP12 in CHO-K1, C6
lysates using MMP12 antibody.