

FYN Rabbit mAb


ABZ2522

Product Overview


Name	FYN Rabbit mAb
Accession(Primary)	P39688
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; FC-1:20-1:50; IP-1:20-1:50

Target information (P39688)

Gene Name	Fyn
About Gene Symbol	



Western blot analysis of Fyn in U87-MG, 3T3 muscle lysates using FYN antibody.



Western blot analysis of Fyn in rat Brain lysates using Fyn antibody.