


p21 Rabbit mAb

ABZ2523

Product Overview

Name	p21 Rabbit mAb
Accession(Primary)	P38936
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000; IF-1:50-1:200

Target information (P38936)

Gene Name CDKN1A**About Gene Symbol**

Western blot analysis of p21 in MCF-7 lysates using p21 antibody.

