

eIF4A3 Rabbit mAb


ABZ2524

Product Overview


Name	eIF4A3 Rabbit mAb
Accession(Primary)	P38919
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; ICC/IF-1:200-1:500; FC-1:100-1:200; IP-1:20-1:50

Target information (P38919)


Gene Name	EIF4A3
About Gene Symbol	



Western blot analysis of eIF4A3 in HeLa, CHO-K1, C6 lysates using eIF4A3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Eif4a3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Eif4a3 (green)
in HeLa using Eif4a3 antibody, and DAPI (blue)