

alpha Synuclein Rabbit mAb


ABZ2528

Product Overview


Name	alpha Synuclein Rabbit mAb
Accession(Primary)	P37840
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:1000-1:5000; IP-1:20-1:50

Target information (P37840)

Gene Name	SNCA
About Gene Symbol	



Western blot analysis of Alphasynuclein in rat Brain lysates using alpha Synuclein antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Alphasynuclein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.