

PGM1 Rabbit mAb


ABZ2531

Product Overview


Name	PGM1 Rabbit mAb
Accession(Primary)	P36871
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (P36871)

Gene Name	PGM1
About Gene Symbol	



Western blot analysis of PGM1 in C6, 3T3, HeLa lysates using PGM1 antibody.



Immunocytochemistry analysis of PGM1(green) in HeLa using PGM1 antibody, and DAPI(blue)