

MEK2 Rabbit mAb


ABZ2533

Product Overview


Name	MEK2 Rabbit mAb
Accession(Primary)	P36507
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (P36507)


Gene Name	MAP2K2
About Gene Symbol	MAP2K2 (Mitogen-Activated Protein Kinase Kinase 2) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of MEK2 in mouse brain, mouse lung lysates using MEK2 antibody.



Western blot analysis of MEK2 in HeLa, A549 lysates using MEK2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using MEK2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of MEK2 (green) in hela using MEK2 antibody, and DAPI (blue).