

Aryl hydrocarbon Receptor Rabbit mAb


ABZ2534

Product Overview


Name	Aryl hydrocarbon Receptor Rabbit mAb
Accession(Primary)	P35869
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 96 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100-1:200

Target information (P35869)

Gene Name	AHR
About Gene Symbol	AHR (Aryl Hydrocarbon Receptor) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Aryl hydrocarbon Receptor in HeLa lysates using Aryl Hydrocarbon Receptor antibody.



Immunocytochemistry analysis of Aryl hydrocarbon Receptor (green) in HeLa using Aryl hydrocarbon Receptor antibody, and DAPI (blue)