


Sarcomeric alpha Actinin Rabbit mAb

ABZ2537


Product Overview

Name	Sarcomeric alpha Actinin Rabbit mAb
Accession(Primary)	P35609
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 104 kDa; Observed MW: 104 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20


Target information (P35609)

Gene Name ACTN2**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded mouse colon using Sarcomeric alpha Actinin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Sarcomeric Alpha Actinin in mouse skeletal muscle, mouse heart lysates using Sarcomeric Alpha Actinin antibody.



Western blot analysis of Sarcomeric Alpha Actinin in rat Brain lysates using Sarcomeric Alpha Actinin antibody