

MYH10 Rabbit mAb


ABZ2538

Product Overview


Name	MYH10 Rabbit mAb
Accession(Primary)	P35580
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 229 kDa; Observed MW: 229 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (P35580)

Gene Name	MYH10
About Gene Symbol	



Western blot analysis of nonmuscle Myosin IIB/MYH10 in 3T3 lysates using MYH10 antibody.



Western blot analysis of nonmuscle Myosin IIB/MYH10 in A549, HL-60, U251, U87-MG lysates using nonmuscle Myosin IIB/MYH10 antibody.