

MYH9 Rabbit mAb


ABZ2539

Product Overview


Name	MYH9 Rabbit mAb
Accession(Primary)	P35579
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 227 kDa; Observed MW: 227 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IP-1:20-1:50

Target information (P35579)

Gene Name	MYH9
About Gene Symbol	MYH9 (Myosin Heavy Chain 9) is a structural cytoskeletal protein involved in adaptive immunity. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of nonmuscle Myosin IIA in HeLa, A549, HL-60, U2OS, C6 lysates using MYH9 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using nonmuscle Myosin IIA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.