

IRS1 Rabbit mAb


ABZ2540

Product Overview

Name	IRS1 Rabbit mAb
Accession(Primary)	P35568
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 132 kDa; Observed MW: 170 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:20-1:50; ICC/IF-1:50-1:200; FC-1:20


Target information (P35568)

Gene Name	IRS1
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of IRS1 in 293 lysates using IRS1 antibody.



Western blot analysis of IRS1 in HEK293 lysates using IRS1 antibody.