

MAS1L Rabbit mAb

ABZ2542


Product Overview

Name	MAS1L Rabbit mAb
Accession(Primary)	P35410
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100; IP-1:100-1:500

Target information (P35410)

Gene Name MAS1L

About Gene Symbol



Western blot analysis of MAS1L in K562, rat Brain, C6, 3T3, HeLa lysates using MAS1L antibody.

