

ICAM2 Rabbit mAb


ABZ2545

Product Overview


Name	ICAM2 Rabbit mAb
Accession(Primary)	P35330
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 31 kDa; Observed MW: 50-60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:1000-1:5000; IP-1:20-1:50

Target information (P35330)

Gene Name	Icam2
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse lung using ICAM2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of ICAM2 in HeLa, A549, HL-60, U2OS lysates using ICAM2 antibody.