

Prohibitin Rabbit mAb



ABZ2548

Product Overview

Name	Prohibitin Rabbit mAb
Accession(Primary)	P35232
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; IP-1:50-1:100


Target information (P35232)

Gene Name	PHB1
About Gene Symbol	PHB1 (Prohibitin 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Prohibitin (green) in HeLa using Prohibitin antibody, and DAPI (blue).

Western blot analysis of Prohibitin in HeLa, A549, HL-60, U2OS, C6 lysates using Prohibitin antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Prohibitin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.