

Bmi1 Rabbit mAb


ABZ2549

Product Overview


Name	Bmi1 Rabbit mAb
Accession(Primary)	P35226
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 37 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:20-1:50; FC-1:50-1:100; IP-1:20-1:50

Target information (P35226)


Gene Name	BMI1
About Gene Symbol	BMI1 (BMI1 Proto-Oncogene, Polycomb Ring Finger) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targe



Immunocytochemistry analysis of BMI1 (green) in Hela using BMI1 antibody, and DAPI (blue).



Western blot analysis of Bmi1 in 3T3 lysates using Bmi1 antibody.



Western blot analysis of Bmi1 in C6, rat Brain, Ramos lysates using Bmi1 antibody