

beta Catenin Rabbit mAb


ABZ2551

Product Overview


Name	beta Catenin Rabbit mAb
Accession(Primary)	P35222
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:1000-1:20000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:20

Target information (P35222)


Gene Name	CTNNB1
About Gene Symbol	



Immunocytochemistry analysis of beta Catenin (green) in 293T using beta Catenin antibody, and DAPI(blue).



Western blot analysis of beta Catenin in NTERA-2, 3T3 lysates using beta Catenin antibody.



Western blot analysis of beta Catenin in C6
lysates using beta Catenin antibody