

# Glypican 1/ GPC1 Rabbit mAb


ABZ2553

## Product Overview


|                             |                                                       |
|-----------------------------|-------------------------------------------------------|
| <b>Name</b>                 | Glypican 1/ GPC1 Rabbit mAb                           |
| <b>Accession(Primary)</b>   | P35052                                                |
| <b>Host Species</b>         | Rabbit                                                |
| <b>Clonality</b>            | Monoclonal                                            |
| <b>Isotype</b>              | IgG                                                   |
| <b>Species Reactivity</b>   | Human                                                 |
| <b>Calculated MW</b>        | Calculated MW: 62 kDa; Observed MW: 55 kDa            |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; ICC/IF-1:200-1:500;<br>FC-1:50-1:100 |

## Target information (P35052)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | GPC1 |
| <b>About Gene Symbol</b> |      |



Western blot analysis of Glypican 1/GPC1 in Jurkat, HeLa lysates using Glypican 1 antibody.



Immunocytochemistry analysis of Glypican 1/GPC1 (green) in HeLa using Glypican 1/GPC1 antibody, and DAPI (blue)