

GTF2H1 (3E8) Mouse mAb

ABZ2554


Product Overview

Name	GTF2H1 (3E8) Mouse mAb
Accession(Primary)	P32780
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 62 kDa; Observed MW: 62 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20


Target information (P32780)

Gene Name GTF2H1


About Gene Symbol



Western blot analysis of GTF2H1 in HeLa, HeLa NE and Jurkat lysates using GTF2H1 antibody.



Immunocytochemistry analysis of GTF2H1 (3E8) in HeLa using GTF2H1 antibody.



Immunoprecipitation analysis of GTF2H1 (3E8) in HeLa lysates using GTF2H1 antibody.