


MCM4 Rabbit mAb

ABZ2557

Product Overview

Name	MCM4 Rabbit mAb
Accession(Primary)	P33991
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 97 kDa; Observed MW: 97 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:500

Target information (P33991)

Gene Name MCM4**About Gene Symbol**

Western blot analysis of MCM4 in HeLa, A549, MCF-7, U251 lysates using MCM4 antibody.

