

kappa Opioid Receptor Rabbit mAb

ABZ2558


Product Overview

Name	kappa Opioid Receptor Rabbit mAb
Accession(Primary)	P33534
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; IP-1:20-1:50

Target information (P33534)

Gene Name Oprk1

About Gene Symbol



Western blot analysis of Kappa Opioid Receptor in rat Brain lysates using kappa Opioid Receptor antibody.

