

GBP1 Rabbit mAb


ABZ2560

Product Overview


Name	GBP1 Rabbit mAb
Accession(Primary)	P32455
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200

Target information (P32455)


Gene Name	GBP1
About Gene Symbol	



Immunocytochemistry analysis of GBP1 (green) in Hela using GBP1 antibody, and DAPI (blue).



Western blot analysis of GBP1 in HepG2, U251, HL-60, Lncap lysates using GBP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using GBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.