

## CXCR5 Rabbit mAb

ABZ2561


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | CXCR5 Rabbit mAb                           |
| <b>Accession(Primary)</b>   | P32302                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 42 kDa; Observed MW: 42 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; ICC/IF-1:50-1:100         |

### Target information (P32302)

**Gene Name** CXCR5

**About Gene Symbol**



Western blot analysis of CXCR5 in K562, C6, 3T3, HeLa lysates using CXCR5 antibody.

