

CCR7 Rabbit mAb


ABZ2562

Product Overview

Name	CCR7 Rabbit mAb
Accession(Primary)	P32248
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; IP-1:20

Target information (P32248)

Gene Name	CCR7
About Gene Symbol	



Western blot analysis of CCR7 in rat Brain, C6, 3T3, Hela lysates using CCR7 antibody.