

Peroxiredoxin 2 Rabbit mAb


ABZ2564

Product Overview


Name	Peroxiredoxin 2 Rabbit mAb
Accession(Primary)	P32119
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:500; FC-1:200-1:1000

Target information (P32119)

Gene Name	PRDX2
About Gene Symbol	



Western blot analysis of Peroxiredoxin 2/PRP in K562 lysates using Peroxiredoxin 2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Peroxiredoxin 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.