

Peroxiredoxin 2 Rabbit mAb


ABZ2565

Product Overview


Name	Peroxiredoxin 2 Rabbit mAb
Accession(Primary)	P32119
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:100-1:1000

Target information (P32119)

Gene Name	PRDX2
About Gene Symbol	



Western blot analysis of Peroxiredoxin 2/PRP in 3T3, mouse brain lysates using Peroxiredoxin 2 antibody.



Western blot analysis of Peroxiredoxin 2/PRP in HeLa, A549, C6 lysates using Peroxiredoxin 2/PRP antibody.