

L1CAM Rabbit mAb


ABZ2566

Product Overview

Name	L1CAM Rabbit mAb
Accession(Primary)	P32004
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 140 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (P32004)

Gene Name	L1CAM
About Gene Symbol	L1CAM (L1 Cell Adhesion Molecule) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of L1CAM in A375 lysates
using L1CAM antibody.