


## CD32B Rabbit mAb

ABZ2567


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | CD32B Rabbit mAb                           |
| <b>Accession(Primary)</b>   | P31994                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 34 kDa; Observed MW: 44 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:1000       |

### Target information (P31994)

**Gene Name** FCGR2B**About Gene Symbol**

Western blot analysis of CD32B in C6 lysates  
using CD32B antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CD32B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.