

STIP1/STI1 Rabbit mAb


ABZ2568

Product Overview


Name	STIP1/STI1 Rabbit mAb
Accession(Primary)	P31948
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 63 kDa; Observed MW: 63 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; FC-1:20-1:50

Target information (P31948)


Gene Name	STIP1
About Gene Symbol	



Western blot analysis of STIP1 in mouse brain, mouse heart lysates using STIP1 antibody.



Western blot analysis of STIP1 in K562, rat Brain, C6 lysates using STIP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using STIP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.