

## 14-3-3 alpha/beta Rabbit mAb


ABZ2569

### Product Overview


|                             |                                                           |
|-----------------------------|-----------------------------------------------------------|
| <b>Name</b>                 | 14-3-3 alpha/beta Rabbit mAb                              |
| <b>Accession(Primary)</b>   | P31946                                                    |
| <b>Host Species</b>         | Rabbit                                                    |
| <b>Clonality</b>            | Monoclonal                                                |
| <b>Isotype</b>              | IgG                                                       |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                         |
| <b>Calculated MW</b>        | Calculated MW: 28 kDa; Observed MW: 28 kDa                |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC:1:200-1:1000;<br>ICC:/IF-1:100-1:200 |

### Target information (P31946)


|                          |                                                                                                                                                                                                                                                                                        |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | YWHAB                                                                                                                                                                                                                                                                                  |
| <b>About Gene Symbol</b> | YWHAB is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of 14-3-3 alpha/beta in rat brain, rat lung lysates using 14-3-3 alpha/beta antibody.



Western blot analysis of 14-3-3 alpha/beta in HeLa, A549, MCF-7, U251, C2C12 lysates using 14-3-3 alpha/beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using 14-3-3 alpha/beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.