


S100A9 Rabbit mAb

ABZ2572

Product Overview

Name	S100A9 Rabbit mAb
Accession(Primary)	P31725
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 13 kDa; Observed MW: 13 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P31725)

Gene Name S100A9**About Gene Symbol**

Western blot analysis of S100A9 in mouse spleen, mouse lung lysates using S100 alpha9 antibody.

