

TIA1 Rabbit mAb


ABZ2573

Product Overview


Name	TIA1 Rabbit mAb
Accession(Primary)	P31483
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; IP-1:10-1:100

Target information (P31483)

Gene Name	TIA1
About Gene Symbol	TIA1 (TIA1 Cytotoxic Granule Associated RNA Binding Protein) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to



Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of TIA1 in (1) Jurkat lysates;
(2) NIH/3T3 lysates using TIA1 antibody.