

RRM2 Rabbit mAb


ABZ2574

Product Overview


Name	RRM2 Rabbit mAb
Accession(Primary)	P31350
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; IP-1:10-1:100

Target information (P31350)

Gene Name	RRM2
About Gene Symbol	



Immunocytochemistry analysis of RRM2 (green) in Jurkat using RRM2 antibody, and DAPI (blue).



Western blot analysis of RRM2 in K562, HeLa lysates using RRM2 antibody.