

SDHA Rabbit mAb


ABZ2576

Product Overview


Name	SDHA Rabbit mAb
Accession(Primary)	P31040
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:500; ICC/IF-1:100-1:200; FC-1:20; IP-1:10-1:50

Target information (P31040)


Gene Name	SDHA
About Gene Symbol	



Western blot analysis of SDHA in K562, C6, 3T3 lysates using SDHA antibody.



Immunocytochemistry analysis of SDHA(green) in Hela using SDHA antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using SDHA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.