

## CLIP170 Rabbit mAb


ABZ2577

### Product Overview


|                             |                                                           |
|-----------------------------|-----------------------------------------------------------|
| <b>Name</b>                 | CLIP170 Rabbit mAb                                        |
| <b>Accession(Primary)</b>   | P30622                                                    |
| <b>Host Species</b>         | Rabbit                                                    |
| <b>Clonality</b>            | Monoclonal                                                |
| <b>Isotype</b>              | IgG                                                       |
| <b>Species Reactivity</b>   | Human                                                     |
| <b>Calculated MW</b>        | Calculated MW: 162 kDa; Observed MW: 170 kDa              |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500;<br>ICC/IF-1:50-1:100 |

### Target information (P30622)


|                          |                                                                                                                                                                                                                                                                                                                          |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CLIP1                                                                                                                                                                                                                                                                                                                    |
| <b>About Gene Symbol</b> | CLIP1 (CAP-Gly Domain Containing Linker Protein 1) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of CLIP170 in HeLa, A549, HL-60, U251, U87-MG lysates using CLIP170 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CLIP170 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of CLIP17 (green)  
in hela using CLIP170 antibody, and DAPI(blue)