


Angiotensin II Type 1 Receptor Rabbit mAb

ABZ2578

Product Overview

Name	Angiotensin II Type 1 Receptor Rabbit mAb
Accession(Primary)	P30556
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P30556)

Gene Name AGTR1**About Gene Symbol**

Western blot analysis of Angiotensin II Type 1 Receptor in HepG2 lysates using Angiotensin II Type 1 Receptor antibody.

